# Emicizumab Calibrator

## EMICIZUMAB PLASMA CALIBRATOR

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### ENGLISH VERSION INTENDED USE

**Emicizumab Calibrator** is a calibrator plasma intended for use in the calibration of a modified one-stage FVIII activity assay for determination of active emicizumab concentration from a clot time measurement performed on an automated coagulation analyzer.

### SUMMARY

Emicizumab (1), a drug for the prophylactic treatment of patients with hemophilia A, is a bispecific antibody that bridges activated factor IX (FIXa) and factor X (FX) to restore the function of missing FVIII that is needed for effective hemostasis. **Emicizumab Calibrator** may be used with a modified one-stage FVIII activity assay on an automated coagulation analyzer to determine the active emicizumab concentration in citrated plasma.

### PRINCIPLE OF THE PROCEDURE

**Emicizumab Calibrator** enables the quantification, of the degree of correction of the APTT based one-stage assay by emicizumab present in the patient’s sample using a coagulation instrument. The calibrator is diluted to generate a series of concentration levels and each level is combined with FVIII deficient plasma, APTT reagent and incubated. The reaction is initiated by addition of CaCl<sub>2</sub> and the clot time is measured. The degree of correction of the APTT is related to the emicizumab activity and is converted into µg/mL using the calibration curve.

### REAGENTS

**Catalog No: 151-201-CE**

**Contents: 5 x 1.0 mL vials**

**Ingredients: Emicizumab Calibrator** is prepared from FVIII immunodepleted citrated human plasma spiked with 100µg/mL emicizumab. The plasma calibrator contains buffer and stabilizers and is lyophilized. **Emicizumab Calibrator** does not contain any preservatives.

**Preparation for Use: Emicizumab Calibrator** should be reconstituted with 1mL of deionized or distilled water. Swirl gently; do not shake. Allow the calibrator to stand for 15 minutes at room temperature to ensure complete dissolution before use.

**Storage and Stability:** Unopened **Emicizumab Calibrator** is stable until the expiration date shown on the label when stored at 2-8°C. Reconstituted control material is stable (e.g., less than a 10% shift in the baseline values) for 8 hours when stored capped at 2-8°C or 4 hours at room temperature (23-25°C).

### WARNINGS

Standard precautions should always be taken, including wearing personal protective equipment, when handling potentially biohazardous human sourced material such as citrated plasma. **Emicizumab Calibrator** material is prepared with plasma collected from donors screened for CJD and which was tested at source and found negative for HBsAG, syphilis and antibodies to HIV and HCV and nonreactive for HIV-1 RNA and HCV RNA by FDA approved tests.

All biohazardous material including **Emicizumab Calibrator** should be disposed of according to current local, state and federal regulations.

### TRACEABILITY

**Emicizumab Calibrator** lots are assigned using an r<sup>2</sup> Diagnostics master calibrator lot that has been assigned based on the manufacturer’s potency assignment for the drug.

### PROCEDURE

Testing is performed on a coagulation analyzer using mechanical or optical based clot detection. Individual system procedural details (e.g. reagents, volume, dilutions, incubation times) may be found in the instrument specific application sheets.

### EXPECTED VALUES

Emicizumab is spiked into FVIII deficient plasma at 100µg/mL and the actual active emicizumab concentration is assigned versus a master calibrator lot. The exact concentration of emicizumab in each calibrator vial is indicated on the vial label and the box’s outer label in each kit.

### LIMITATIONS

Performance of **Emicizumab Calibrator** should be verified under individual laboratory conditions using **Emicizumab Controls Level 1 and Level 2**.

FVIII concentrate replacement products present in patient samples may cause an elevation in emicizumab concentrations using the modified FVIII one-stage assay.

### REFERENCES

- ↑ Kitazawa T, et al, (2012). “A Bispecific Antibody to Factors IXa and X Restores Factor VIII Hemostatic Activity in a Hemophilia A Model”. Nature Medicine 18: 1570-1574.

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