NoFact IX

Immunodepleted Factor IX Deficient Substrate Plasma

RESULTS

Results of a factor assay may be expressed in % activity or IU/ml. The analytical measurement range (linearity) is 1% - 188% FIX activity.

LIMITATIONS

Hemolysis to 500 mg/dL hemoglobin, icterus to 20 mg/dL unconjugated bilirubin, and lipemia to 2000 mg/dL triglycerides cause less than a 10% shift in % FIX recoveries using NoFACT IX Deficient Plasma with Stago PTT-A on the Stago Compact. Unfractionated heparin, Low Molecular Weight Heparin, and direct thrombin inhibitors interfere with FIX determinations. Lupus anticoagulants may also interfere (7).

The performance characteristics of NoFACT IX Deficient Plasma were not evaluated for other coagulation analyzers and APTT reagents or coagulation systems.

PERFORMANCE CHARACTERISTICS

The method comparison and analytical studies of NoFACT IX Deficient Plasma were assessed using Diagnostica Stago STA Compact coagulation analyzers, Stago PTT Automate 5 (PTT-A), and Stago Unicibrator and controls.

Precision:

CLSI EP5-A2 (5) precision estimates of the Stago PTT-A Factor IX assay using NoFACT IX Deficient Plasma, as % CV of the recovered FIX values, were:

Plasma	Mean FIX activity, %	% CV, Within-run (S-r)	% CV, Lot-to-Lot (S-lot)	% CV, Within-Device (S-device)
System N (Normal Control Plasma) N=240	109%	5.1%	1.1%	7.2%
System P (Abnormal Control Plasma) N=240	43%	4.6%	2.9%	7.7%
Low FIX pooled patient plasma N=120	14%	5.7%	2.2%	7.8%

Correlation:

Two hundred and thirty-three plasma samples from patients and donors were assessed in three laboratories in parallel with the Stago PTT-A FIX assay using Stago IX Deficient plasma and with the Stago PTT-A FIX assay using NoFACT IX Deficient Plasma. The regression statistics were:

	All Labs N=233	Site 1 N=100	Site 2 N=80	Site 3 N=53
Pente	0.858	0.785	0.923	0.817
Ordonnée	5.729	6.168	5.328	11.206
r ²	0.915	0.977	0.907	0.830
r	0.956	0.988	0.988	0.911

REFERENCE VALUES

A typical reference range is 78%-184% (6), however each laboratory should determine a reference range for FIX activity for its particular population and instrument / reagent system.

LITERATURE REFERENCES

- Bajaj S and Thompson A, "Molecular and Structural Biology of Factor IX", in Hemostasis and Thrombosis: Basic Principles and Clinical Practice, 5th edition, editors Coleman R et al, Lippincott Williams and Wilkins, 2006.
- Marques MB and Fritsma GA, Quick Guide to Coagulation Testing, AAC Press2006.
- Laposata M, et al., The Clinical Hemostasis Handbook, Year Book Medical Publishers, 1989.
- H21-A5, Collection, Transport and Processing of Blood Samples for Testing Plasma-based Coagulation Assays and Molecular Hemostasis Assays; Approved Guideline-Fifth Edition, Clinical Laboratory Standards Institute, 2008.
- EP5-A2, Evaluation of Precision Performance of Quantitative Measurement Methods; Approved Guideline-2nd Edition, Clinical Laboratory Standards Institute, 2004.
- Bennett ST, Lehman CM, and Rodgers GM, Laboratory Hemostasis: A Practical Guide for Pathologists, Springer, 2007.
- Brandt JT, et al, Effect of lupus anticoagulants on the activated partial thromboplastin time, Arch Pathol Lab Med 115: 109, 1991.

TEST PROCEDURE

Contact r2 Diagnostics for instrument applications using APTT reagent to test for FIX concentration.

Quality Control

Quality control of coagulation tests involves multiple components. Each laboratory should establish a quality control program that includes both normal and abnormal control plasmas.

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- 5) EP5-A2, Evaluation of Precision Performance of Quantitative Measurement Methods; Approved Guideline-2nd Edition, Clinical Laboratory Standards Institute, 2004.
- 6) Bennett ST, Lehman CM, and Rodgers GM, Laboratory Hemostasis: A Practical Guide for Pathologists, Springer, 2007.
- 7) Brandt JT, et al. Effect of lupus anticoagulants on the activated partial thromboplastin time, Arch Pathol Lab Med 115: 109, 1991

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BIBLIOGRAFIA

- 1) Bajaj S and Thompson A, "Molecular and Structural Biology of Factor IX", in Hemostasis and Thrombosis: Basic Principles and Clinical Practice, 5th edition, editors Coleman R et al, Lippincott Williams and Wilkins, 2006.
- 2) Marques MB and Fritsma GA, Quick Guide to Coagulation Testing, AACC Press2006.
- 3) Laposata M, et al., The Clinical Hemostasis Handbook, Year Book Medical Publishers, 1989.
- 4) H21-A5, Collection, Transport and Processing of Blood Samples for Testing Plasma-based Coagulation Assays and Molecular Hemostasis Assays; Approved Guideline-Fifth Edition, Clinical Laboratory Standards Institute, 2008.
- 5) EP5-A2, Evaluation of Precision Performance of Quantitative Measurement Methods; Approved Guideline-2nd Edition, Clinical Laboratory Standards Institute, 2004.
- 6) Bennett ST, Lehman CM, and Rodgers GM, Laboratory Hemostasis: A Practical Guide for Pathologists, Springer, 2007.
- 7) Brandt JT, et al. Effect of lupus anticoagulants on the activated partial thromboplastin time, Arch Pathol Lab Med 115: 109, 1991

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5)

EP5-A2, Evaluation of Precision Performance of Quantitative Measurement Methods; Approved Guideline-2nd Edition, Clinical Laboratory Standards Institute, 2004.

Bennett ST, Lehman CM, and Rodgers GM, Laboratory Hemostasis: A Practical Guide for Pathologists, Springer, 2007.

6)

Brandt JT, et al. Effect of lupus anticoagulants on the activated partial thromboplastin time, Arch Pathol Lab Med 115: 109, 1991

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- 1) Bajaj S and Thompson A, "Molecular and Structural Biology of Factor IX", in Hemostasis and Thrombosis: Basic Principles and Clinical Practice, 5th edition, editors Coleman R et al, Lippincott Williams and Wilkins, 2006.
- 2) Marques MB and Fritsma GA, Quick Guide to Coagulation Testing, AACC Press2006.
- 3) Laposata M, et al., The Clinical Hemostasis Handbook, Year Book Medical Publishers, 1989.
- 4) H21-A5, Collection, Transport and Processing of Blood Samples for Testing Plasma-based Coagulation Assays and Molecular Hemostasis Assays; Approved Guideline-Fifth Edition, Clinical Laboratory Standards Institute, 2008.