PHOSPHOPLASTIN RL

Prothrombin Time Reagent

VERSION FRANCAISE
APPLICATION

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LITERATUR

- Quick A.J., The Prothrombin Time in Hemophilia and in Obstructive Jaundice. J. Biol. Chem.: 109,73-74; 1935.
- Biggs R. ed., Human Blood Coagulation Hemostasis and Thrombosis Second Ed. Blackwell Scientific Publications, London 1976.
- Peterson C.E., Kwaan H.C., Current Concepts of Warfarin Therapy, Arch Intern. Med.146: 581-584, 1986.
- Loeliger E.A.: ICEH/ICTH Recommendations for Reporting Prothrombin Time in Oral Anticoagulant Control, Throm. Haemost. 53: 155-156, 1985.

ENGLISH VERSION

INTENDED USE

Phosphoplastin RL is intended for use in the one-stage prothrombin time (PT) test and in PT based hemostasis assays. The PT can be used to screen for coagulopathies, for extrinsic factor assays, and for therapeutic monitoring of Warfarin and related drugs.

PRINCIPLE

Tissue thromboplastin in the presence of calcium and FVIIa activates the tissue factor (extrinsic) pathway of coagulation. Deficiencies in the extrinsic and common pathways prolong the PT.

COMPOSITION

Phosphoplastin RL consists of a liquid saline extract of rabbit brain, calcium ions, preservatives and stabilizers.

PREPARATION FOR USE

Phosphoplastin RL comes ready for use. The vial should be inverted several times until a homogeneous suspension is obtained.

STORAGE AND STABILITY

Phosphoplastin RL should be stored at 2 to 8°C and is stable until the expiration date indicated on the vial. After opening, the original activity is stable for fourteen (14) days when stored at 2 to 8°C. DO NOT FREEZE. The reagent is a fine suspension of rabbit brain particles. Large flaky particles in the suspension or prolonged PT times on testing normal plasma or controls may indicate product deterioration.

SPECIMEN COLLECTION AND HANDLING

Plasma is obtained from whole blood anti-coagulated with 1 part 3.2% sodium citrate to 9 parts whole blood. Process the collected whole blood and handle the plasma according to CLSI guideline H21-A4 (or superseding edition).

MATERIALS AND EQUIPMENT REQUIRED, BUT NOT PROVIDED

- Coagulation Instrument or 37°C Water Bath and timer
- Reaction Cups or Plastic test tubes (12 x 75 mm)
- Pipetting device to deliver 200 µL or 100 µL.
- Centrifuge
- Distilled or de-ionized water
- Control plasmas such as r2 Diagnostics PlasmaCon N (30-201), L1 (31-201), and L2 (32-201).

TEST PROCEDURE FOR PT

NOTE: Throughout the procedure, all test tubes, syringes, and pipettes should be plastic.

I. Automated And Semi-Automated Methods

Refer to the appropriate Instrument Operator’s Manual for detailed instructions.

II. Manual Method

- Collect and prepare the blood specimen according to the directions outlined in SPECIMEN COLLECTION AND HANDLING.
- Reconstitute the control plasmas according to their package inserts.
- Perform all tests in duplicate.
- Prewarm the PT reagent to 37°C for at least 10 minutes.
- Prewarm 100 µL of the test plasma or control plasma for 2-3 minutes at 37°C.
- Add 200µL **Phosphoplastin RL** reagent to the plasma, simultaneously starting a stopwatch and record the time required for clot formation in seconds.

QUALITY CONTROL

Each laboratory should establish a quality control program that includes normal and abnormal controls to evaluate instrument, reagent and technologist performance. The normal and abnormal controls should be tested daily prior to performing tests on patient plasmas. Control ranges should be reestablished upon any change in lot number of reagent or control, or upon major instrument service. If the controls do not recover within their reference ranges, patient results should be considered invalid and not reported.

RESULTS

The results of the Prothrombin Time tests can be reported in units of seconds, % activity, or INR. Results greater than the upper limits of the range should be considered abnormal and follow-up testing should be performed. PT values below the lower limits of the range may indicate a compromised sample

or instrument and should be investigated.

LIMITATIONS

Phosphoplastin RL is for *in vitro* diagnostic use only. As

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- Quick A.J., The Prothrombin Time in Hemophilia and in Obstructive Jaundice. J. Biol. Chem.: 109,73-74; 1935.
- Biggs R. ed , Human Blood Coagulation Hemostasis and Thrombosis Second Ed. Blackwell Scientific Publications, London 1976.
- Peterson C.E., Kwaan H.C., Current Concepts of Warfarin Therapy. Arch Intern. Med.146: 581-584, 1986.
- Loeliger E.A.: ICEH/ICTH Recommendations for Reporting Prothrombin Time in Oral Anticoagulant Control, Throm. Haemost. 53: 155-156, 1985.

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© 2011 R2 Diagnostics, Inc.