IVD

PHOSPHOPLASTIN RL

CE

Prothrombin Time Reagent

ENGLISH VERSION

INTENDED USE

Phosphoplastin RL is intended for use in the one-stage prothrombin time (PT) test and in PT based hemostasis assays. The PT can be used to screen for coagulopathies, for extrinsic factor assays, and for therapeutic monitoring of Warfarin and related drugs.

PRINCIPLE

Tissue thromboplastin in the presence of calcium and FVIIa activates the tissue factor (extrinsic) pathway of coagulation. Deficiencies in the extrinsic and common pathways prolong the PT.

COMPOSITION

Phosphoplastin RL consists of a liquid saline extract of rabbit brain, calcium ions, preservatives and stabilizers.

PREPARATION FOR USE

Phosphoplastin RL comes ready for use. The vial should be inverted several times until a homogeneous suspension is obtained.

STORAGE AND STABILITY

Phosphoplastin RL should be stored at 2 to 8°C and is stable until the expiration date indicated on the vial. After opening, the original activity is stable for fourteen (14) days when stored at 2 to 8°C. DO NOT FREEZE. The reagent is a fine suspension of rabbit brain particles. Large flaky particles in the suspension or prolonged PT times on testing normal plasma or controls may indicate product deterioration.

SPECIMEN COLLECTION AND HANDLING

Plasma is obtained from whole blood anti-coagulated with 1 part 3.2% sodium citrate to 9 parts whole blood. Process the collected whole blood and handle the plasma according to CLSI guideline H21-A4 (or superseding edition).

MATERIALS AND EQUIPMENT REQUIRED, BUT NOT PROVIDED

- Coagulation Instrument or 37°C Water Bath and timer
- Reaction Cups or Plastic test tubes (12 x 75 mm)
- Pipetting device to deliver 200 µL or 100 µL.
- Centrifuge
- Distilled or de-ionized water
- Control plasmas such as r2 Diagnostics PlasmaCon N (30-201), L1 (31-201), and L2 (32-201).

TEST PROCEDURE FOR PT

NOTE: Throughout the procedure, all test tubes, syringes, and pipettes should be plastic.

I. Automated And Semi-Automated Methods

Refer to the appropriate Instrument Operator's Manual for detailed instructions.

II. Manual Method

- Collect and prepare the blood specimen according to the directions outlined in SPECIMEN COLLECTION AND HANDLING.
- Reconstitute the control plasmas according to their package inserts.
- Perform all tests in duplicate.
- Prewarm the PT reagent to 37°C for at least 10 minutes.
- Prewarm 100 µL of the test plasma or control plasma for 2-3 minutes at 37°C.
- Add 200µL **Phosphoplastin RL** reagent to the plasma, simultaneously starting a stopwatch and record the time required for clot formation in seconds.

QUALITY CONTROL

Each laboratory should establish a quality control program that includes normal and abnormal controls to evaluate instrument, reagent and technologist performance. The normal and abnormal controls should be tested daily prior to performing tests on patient plasmas. Control ranges should be reestablished upon any change in lot number of reagent or control, or upon major instrument service. If the controls do not recover within their reference ranges, patient results should be considered invalid and not reported.

RESULTS

The results of the Prothrombin Time tests can be reported in units of seconds, % activity, or INR. Results greater than the upper limits of the range should be considered abnormal and follow-up testing should be performed. PT values below the lower limits of the range may indicate a compromised sample

or instrument and should be investigated.

LIMITATIONS

Phosphoplastin RL is for *in vitro* diagnostic use only. As

expected values for the prothrombin time test will vary from one laboratory to another, depending on numerous variables such as the choice of clot detection method or environmental conditions, laboratories should establish their own expected values for patients and well defined ranges for control plasmas. The use of Icteric, Lipemic, or Hemolyzed sample should be avoided due to possible interference especially when using photo-optical instruments. The impact of other therapeutic drugs, in addition to oral anticoagulant therapy, can influence interpretation of PT test results.

EXPECTED VALUES

A reference range study was conducted using triplicate specimens from 120 normal health adults. Approximately equal numbers of males and females were used. The PT results were as follows:

	Mean (seconds)	
	Range for +/- 2 SD	
Photo-optical	12.3	9.7-14.9 seconds
Mechanical	12.6	11.1-14.1 seconds

These values should only serve as guidelines. Because differences may exist between instruments, laboratories, and local populations, it is necessary that each laboratory establish its own reference range of expected prothrombin time results.

DETERMINATION OF INR

The INR is calculated by using the ratio of the Patient PT to the mean of the normal reference range raised to the power of the reagent International Sensitivity Index (ISI).

INR = (Patient PT / mean of normal range)^{ISI}

The ISI is assigned by comparison to a highly sensitive WHO thromboplastin standard reference material. The lower the ISI for a given reagent the more sensitive the reagent is to coagulation factors. The lot specific ISI value for **Phosphoplastin RL** can be found on the outer box front panel.

PERFORMANCE CHARACTERISTICS

I. Precision

Precision studies were performed to establish Within-Run and Between-Run CV's for normal control plasma and abnormal control plasma. A single lot number of **Phosphoplastin RL** reagent was used for these studies. Results are shown below.

Photo-optical clot detection

Control Plasma	Mean PT (sec)	Within Run CV (%)	Between Run CV (%)
PlasmaCon N	12.7	0.7	1.03
PlasmaCon L2	52.9	1.9	1.31

Mechanical clot detection

Control Plasma	Mean PT (sec)	Within Run CV (%)	Between Run CV (%)
PlasmaCon N	14.4	2.1	1.04
PlasmaCon L2	59.7	1.5	0.54

II. Correlation

A correlation study was performed using the **Phosphoplastin RL** Reagent (Y) and another thromboplastin reagent (X). One hundred five (105) plasma specimens, both normal and abnormal clinical samples were tested with both PT reagents. The linear regression equation and coefficient of determination (r²) of the INR values were:

Method	N	R2	Equation
Photo-optical	105	0.98	Y = 1.065X – 0.1245
Mechanical	105	0.92	Y = 1.059X – 0.114

REFERENCE

- Quick A.J., The Prothrombin Time in Hemophilia and in Obstructive Jaundice. J. Biol. Chem.: 109,73-74; 1935.
- Biggs R. ed., Human Blood Coagulation Hemostasis and Thrombosis Second Ed. Blackwell Scientific Publications, London 1976.
- Peterson C.E., Kwaan H.C., Current Concepts of Warfarin Therapy, Arch Intern. Med.146: 581-584, 1986.
- Loeliger E.A.: ICEH/ICTH Recommendations for Reporting Prothrombin Time in Oral Anticoagulant Control, Throm. Haemost. 53: 155-156, 1985.

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1. Quick A.J., The Prothrombin Time in Hemophilia and in Obstructive Jaundice. J. Biol. Chem.: 109,73-74; 1935.
2. Biggs R. ed , Human Blood Coagulation Hemostasis and Thrombosis Second Ed. Blackwell Scientific Publications, London 1976.
3. Peterson C.E., Kwaan H.C., Current Concepts of Warfarin Therapy, Arch Intern. Med.146: 581-584, 1986.
4. Loeliger E.A.: ICEH/ICTH Recommendations for Reporting Prothrombin Time in Oral Anticoagulant Control, Throm. Haemost. 53: 155-156, 1985.

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1. Quick A.J., The Prothrombin Time in Hemophilia and in

R2 Diagnostics, Inc.
1801 Commerce Drive
South Bend, Indiana USA 46628
(574) 288-4377
www.r2diagnostics.com

Obstructive Jaundice. J. Biol. Chem.: 109,73-74; 1935.

2. Biggs R. ed , Human Blood Coagulation Hemostasis and Thrombosis Second Ed. Blackwell Scientific Publications, London 1976.
3. Peterson C.E., Kwaan H.C., Current Concepts of Warfarin Therapy, Arch Intern. Med.146: 581-584, 1986.
4. Loeliger E.A.: ICEH/ICTH Recommendations for Reporting Prothrombin Time in Oral Anticoagulant Control, Throm. Haemost. 53: 155-156, 1985.

Inform of any suspected serious incident related to the device to Info@r2diagnostics.com, or the local r2 distributor, and the competent authority, ministry of health, or delegated agency in the country where the suspected serious incident occurred.

LL-4500 Rev. H
Last Revision: 04/14/25

EU REP

MT Promedt Consulting GmbH
Ernst-Heckel-Straße 7, 66386
St. Ingbert, Germany